

Produced by  prospera

M i n i P e n g u i n M a s s a g e r

P L 0 0 2



O p e r a t i n g M a n u a l

Introduction

Mini Penguin Massager PL002, an easy to use product, delivers a steady relaxing massage to tired muscles. It can relax muscles on your back, shoulder, neck, arms, legs and abdomen with its powerful vibration massage. You can adjust the massage intensity of PL002, which includes two speeds. It is a perfect personal massager for use at home and in office.

Features of PL002 Mini Penguin Massager:

1. Design of ultra-light weight and small size
2. Portable, flexible movement for maximum sensation
3. Muscle relaxation and pain relief on various body areas
4. 1 Button, 2 Vibration speeds

CONTENTS

Safety Warnings... .. 3

Set up... .. 4

Operating the product... .. 5

Trouble shooting... .. 6

Consumer warranty... .. 7

Positions... .. 8

Warranty card... .. 9

Safety Warnings

Consult your doctor prior to using the product if you have heart, skin and/or other health-related issues.

1. This product is intended for adult use only.
2. Do not let babies and children play with the product.
3. Start massaging from the lowest intensity of product. If you feel uncomfortable, reduce the intensity or stop using the product immediately.
4. Do not disassemble the product unless authorized.
5. Do not use the product on eyes and/or any sensitive areas.
6. Do not use the product when the skin is swelling, inflamed, and/or abnormal.

Set up



Unpack the product, and place two AA batteries into the battery compartment on the back of the product.

Product Specifications:

Unit size	15.5 X 5.1 X 4 cm
Unit Weight	84 g (without battery)
Length of the unit	15.5 cm
Supply voltage	Two AA batteries (3V)
Vibration speed	Low : 3600 rpm ; High: 4700 rpm

Operating the product

1. Press the LOW /HIGH /OFF button to turn on the massager, which will start at the low speed of vibration massage.
2. Press the LOW /HIGH /OFF buttons again to adjust the vibration to the higher speed.
3. After the speed is selected, place the massager on the body areas where massage is needed.
4. After the massage is complete, press the LOW /HIGH /OFF button to turn off the massager.

Troubleshooting

Mini Penguin Massager PL002 is a user-friendly product. Follow the operating manual to have it work properly. However, if the following situations happen, please take the recommended troubleshooting procedures:

1. If the vibration head of the product becomes loose, stop using it immediately.

Report the problem to the dealer or Prospera if the product is under warranty.

2. If the product is not working, please check the batteries.

C o n s u m e r W a r r a n t y

Prospera offers a 90-day standard warranty to the massage products. Labor and parts are included during warranty period. Consumers are required to provide model number and dated proof of purchase (sales receipt) when contacting Prospera for the warranty service.

Non-warranty service is provided on a “per incident” basis. The consumer will verify that the product has failed and provide information for servicing the unit. All applicable repairs, parts, shipping, handling, local tax and “per incident” fees will be charged for non-warranty repairs and support calls.

Positions



Warranty Card



Warranty Registration Form

Product Name: _____

Product Model: _____ Purchase Price: _____

Date of Purchase: _____ Retailer: _____

First Name: _____ Last Name: _____

Address: _____

City: _____ State: _____ Zip: _____

Phone: _____ Email: _____

Signature and date: _____

Please mail to:

Prospera Corporation
448 Commerce Way
Livermore, CA 94551

Telephone: 925-225-0888
www.prosperacorp.com
Email:sales@prosperacorp.com

