

Introduction

Prospera's PL016 Personal Jade Massage Cushion is a user friendly product. Prospera's Personal Jade Massage Cushion brings the ultimate power of jade with two pairs of kneading heads and two jade stones. The Personal Jade Massage Cushion can be used on your neck, shoulders, back, thighs, calves and more. Jade is rich in calcium, magnesium and iron—essential minerals for your body. Chinese medicine has long used jade for its healing powers. The kneading mechanism delivers rotating and spot shiatsu massage directly to your muscle tension spots — for instant relief. Each of the four massage heads gives deep tissue Shiatsu massage, helps loosens tight muscles by encouraging blood flow. You can use one of two automated massage programs or apply the massage manually. Speed is adjustable. A perfect personal massager for home and office.

Your muscles will sigh with relief when you soothe them with the Personal Jade Massage Cushion. The ultra-light unit comes with a remote, to bring maximum sensation just where you need it. Firm up muscles at home or in the office!

Features of PL016 Personal Jade Massage Cushion:

1. Multiple speeds to meet your needs
2. Two operating modes: forward and reversal
3. 2 pieces jade stones
4. Deep Shiatsu massage
5. Well designed controller for easy and simple operation
6. A perfect personal massager for the home and office

Product Specifications

- Powerful high-grade motor
- Manual speed regulation
- Two mode selections
- Remote control
- Plugs into power outlet
- Weight 3.5 lbs
- Measures 17.7"x4.5"x13"
- UPC 185277000028

About Prospera

Prospera develops high quality, innovative, affordable massage products — from electronic pulse massagers to jade massagers — for active people. Prospera products are available through leading department stores, consumer electronics stores, and major online retailers. Learn more about Prospera at www.prosperacorp.com

Content

Safety Warnings

Set up

Operating the product

Trouble shooting

Consumer warranty

Warranty card

Prospera company contact

Use positions diagrams

Safety Warnings

When the Personal Jade Massage Cushion is not in use, always unplug it from the transformer. This action conserves energy, prevents electric shock and damage to the product.

Warnings:

Please take these precautions to prevent electric shock, burns, fires, and other injuries.

1. Do not leave this product unattended when it is in use. Unplug the product from the power surge protector or cigarette lighter after usage.
2. Do not operate this product in an excessively humid environment, and avoid direct sunshine upon the product where it is set up.
3. Do not insert any sharp objects into the fabric of the product. These objects might include keys, pens, pencils, kitchen utensils, blades, knives, etc.
4. Do not use the product outdoors.
5. Do not use air refresher spray where the product is set up or in use.
6. Consult your doctor prior to using the product if you have heart, skin and other health issues.
7. This product is not intended for children to use. Keep close supervision if your child uses it.
8. This transformer is recommended to be plugged into a surge protector to avoid damage to its mechanism and electric components.

Attention:

1. Do not let babies and young children to play with the product, remote, or cord.
2. If you feel that the massage is too strong, stop program immediately, and disconnect the product from the transformer.
3. Average daily usage for massage product is about 15 to 30 minutes. Do not use it while you are asleep.
4. This product contains electric components, wires, and other elements that cannot contact excessive moisture. For instance, the bathroom, shower room, steam room, and sauna are places where this product should not be used.
5. Do not sit on the massager to massage.

Set up

Find a good location in your living room , family room , study , or office to set up the Prospera's PL016 Personal Jade Massage Cushion . Use a power surge protector for connecting the product .



Take the product out of retail box , place it at your lower back , or shoulders , and insert the plug to the surge protector . Turn the surge protector on .

Operating the product

Function :

1. Multiple speeds for you to choose
2. Two operating modes: forward and reversal
3. Deep Shiatsu massage

Parameters:

Model: PL016

Name: Personal Jade Massage Cushion

Input: 120V A C , 60H z

Power: < 35 W

Controller



1. Rotation speed is adjustable, low to high
2. Two massage modes: forward and reversal
3. Stop: power off
4. Power indicator is on when massage is in session

Operating instruction:

Turn on the power by choosing massage mode "Forward" or "reversal". Manually regulate speed, turn the dial clockwise to increase the speed for your massage session.

After you are done massaging, switch the Power off to middle to Stop. Turn off the surge protector.

Recommended time of usage: 15 to 30 minutes daily

Trouble shooting

Prospera's PL016 Personal Jade Massage Cushion is a user friendly product. It is built to last. If you follow the operating manual, do not over use it, and follow the precautions, this product will last for a long time.

Should the following situations happen, please take recommended trouble shooting procedures:

1. If the fabric is torn, stop using it immediately.
2. If a wire is broken or the product does not work, unplug the product from the surge protector. Contact Prospera if it is under warranty.

Consumer Warranty

Prospera offers a 90-day standard warranty to the massage products. Labor and parts are included during warranty period. Consumers are required to provide model number and dated proof of purchase (sales receipt) when contacting Prospera for the warranty service.

Non-warranty service is provided on a "per incident" basis. The consumer will verify that the product has failed and provide information for servicing the unit. All applicable repairs, parts, shipping, handling, local tax and "per incident" fees will be charged for non-warranty repairs and support calls.

Warranty card

Prospera and logo

Warranty Registration Form

Product Name: _____

Product Model: _____ Purchase Price: _____

Date of Purchase: _____

Retailer: _____

First Name: _____ Last Name: _____

Address: _____

City: _____ State: _____ Zip: _____

Phone: _____

Email: _____

Signature and date: _____

Please mail to:

Prospera Corporation
405 Boulder Ct, Ste 500
Pleasanton, CA 94566

Telephone: 925-225-0888

Fax: 925-225-0660

Email: sales@prosperacorp.com

Positions

Please include the 4 images Ellen sent to you for positions